

# Estimating $\delta^{13}\text{C}$ of grassland soil-respired $\text{CO}_2$ using an optimized automated chamber system and novel $\text{CO}_2$ isotope Analyzer


Graham Leggett, Richard Vath, Jason Hupp, Doug Lynch, Alex Frodyma, Mark Johnson, Israel Begashaw  
LI-COR Biosciences, Lincoln, Nebraska, USA

## Introduction


- Soil carbon dynamics are affected by various biological and geochemical processes and tracking the evolution of soil-respired  $\text{CO}_2$  provides a powerful solution toward understanding carbon cycling in both natural and managed ecosystems.
- The stable carbon isotope ratio ( $\delta^{13}\text{C}$ ) of soil-respired  $\text{CO}_2$  is an isotopic “fingerprint” which is dynamically affected by soil characteristics, above- and below-ground biological activity, and climate variability.
- This study demonstrates application of a novel, automated chamber-based soil gas measurement system to monitor  $\delta^{13}\text{C}$  of soil-respired  $\text{CO}_2$  in a managed grassland. The Keeling approach was used to estimate source  $\delta^{13}\text{C}$  of respired  $\text{CO}_2$  from natural and manipulated plots, revealing variation due to substrate source and diurnal cycling.

## Methods

Four automated soil gas flux chambers were installed in a managed grassland site near Lincoln, NE, USA in early June 2024. On July 3, 2024, two of the chambers were moved on to sealed collars, one containing *Glycine max* (soybean,  $\text{C}_3$  photosynthetic physiology) residue and one containing *Zea mays* (maize,  $\text{C}_4$ ).




**Figure 1:** Schematic of automated soil gas measurement system. Soil gas flux chambers are connected to a Multiplexer which manages gas sampling and data logging. An LI-7825  $\text{CO}_2$  Isotope /  $\text{NH}_3$  Trace Gas Analyzer interfaces with the multiplexer via a calibration manifold. The manifold provides a novel, passive  $\text{H}_2\text{O}$  equilibration system and facilitates automated daily baseline calibration. Other trace gas analyzers can be connected as well to simultaneously measure fluxes of gases such as  $\text{N}_2\text{O}$  and  $\text{CH}_4$ .




**Figure 2:** Photograph of measurement system

## Results




**Figure 3:** Measured  $\delta^{13}\text{C}$  values for four automated chambers. Week 1 shows data collected over untreated grassland plots while week 2 shows one plot dosed with  $\text{C}_3$  soybean residue (blue) and another with  $\text{C}_4$  maize residue (green). Grey points remained on grassland. The observed  $\delta^{13}\text{C}$  values are consistent with those associated with  $\text{C}_3$  and  $\text{C}_4$  plants: approximately  $-28\text{‰}$  and  $-13\text{‰}$  respectively.

## Estimation of Soil $\delta^{13}\text{C}$



Download  
SoilFluxPro



**Figure 6:** Estimated  $\delta^{13}\text{C}$  from sealed volumes containing soybean (blue) and maize (green) residue, with no spectroscopic correction for  $\text{NH}_3$ . The  $\text{C}_4$  maize off-gassed  $\text{NH}_3$  during decomposition, skewing  $\delta^{13}\text{C}$  to unrealistic values. The LI-7825 analyzer used in this study quantifies  $\text{NH}_3$  to report  $\text{CO}_2$  isotopologues with high accuracy.

## Conclusions

- An automated soil gas measurement system facilitates measurement of  $\text{CO}_2$  fluxes and  $\delta^{13}\text{C}$  over long time periods in field conditions, allowing investigation of changes in carbon cycling due to changes in biological processes.
- A calibration manifold containing a novel, passive water vapor transient equilibration system and automated analyzer baseline system enables long-term, low-power, unattended soil  $\delta^{13}\text{C}$  measurements with minimal maintenance and consumables.
- A new trace gas analyzer provides high-precision measurements of  $\text{CO}_2$  isotopologues in field conditions.